

---

**Product Name: Hexokinase II Rabbit Monoclonal Antibody****Catalog #: AMRe02080**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC,IP
<b>Reactivity</b>	Human,Rat,Hamster
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	0.53mg/ml. The concentration of this product may be batch-dependent.
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
<b>Molecular Weight</b>	Calculated MW: 102 kDa; Observed MW: 102 kDa

**Antigen Information**

<b>Gene Name</b>	HK2
<b>Alternative Names</b>	HKII; HXK2
<b>Gene ID</b>	3099
<b>SwissProt ID</b>	P52789
<b>Immunogen</b>	Recombinant protein of human Hexokinase II


**Background**

MiscellaneousIn vertebrates there are four major glucose-phosphorylating isoenzymes, designated hexokinase I, II, III and IV (glucokinase).


## Research Area

Signal Transduction


## Image Data



Western blot analysis of Hexokinase II in rat skeletal muscle lysates using Hexokinase II antibody.



Western blot analysis of Hexokinase II in HeLa, CHO-K1 lysates using Hexokinase II antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Hexokinase II antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.