

Product Name: HLA C Rabbit Monoclonal Antibody**Catalog #: AMRe02095**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:20-1:50
Molecular Weight	Calculated MW: 41 kDa; Observed MW: 41 kDa

Antigen Information

Gene Name	HLA-C
Alternative Names	MHC; HLAC; HLC-C; D6S204; PSORS1; HLA-JY3
Gene ID	3107.0
SwissProt ID	P10321
Immunogen	


Background

Antigen-presenting major histocompatibility complex class I (MHC I) molecule with an important role in reproduction and antiviral immunity.

Research Area

Immunology

Image Data



Western blot analysis of HLAC in HeLa, 293T, Jurkat, 3T3, HepG2 lysates using HLA C antibody.