

Product Name: HMGCR Rabbit Monoclonal antibody
Catalog #: AMRe02100



Summary

Production Name	HMGCR Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	HMGCR
Alternative Names	LDLCQ3
Gene ID	3156
SwissProt ID	P04035.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 97 kDa; Observed MW: 97 kDa

Background

Transmembrane glycoprotein that is the rate-limiting enzyme in cholesterol biosynthesis as well as in the biosynthesis of

Product Name: HMGCR Rabbit Monoclonal antibody
Catalog #: AMRe02100




nonsterol isoprenoids that are essential for normal cell function including ubiquinone and geranylgeranyl proteins.

Research Area

Cardiovascular

Image Data



Western blot analysis of HMGCR in K562, 3T3, HeLa lysates using HMGCR antibody.

Note

For research use only.