
Product Name: hnRNP K Rabbit Monoclonal Antibody**Catalog #: AMRe02106**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.64mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100
Molecular Weight	Calculated MW: 51 kDa; Observed MW: 62 kDa

Antigen Information

Gene Name	HNRNPK
Alternative Names	HNRNPK; HNRPK; Heterogeneous nuclear ribonucleoprotein K; hnRNP K; Transformation up-regulated nuclear protein; TUNP
Gene ID	3190
SwissProt ID	P61978
Immunogen	A synthetic peptide of human hnRNP K

Background


Facilitate pre-mRNA processing and transport of mRNA from the nucleus to cytoplasm. hnRNP K contains three unique

structural motifs termed KH domains that bind poly(C) DNA and RNA sequences. Intricate architecture enables hnRNP K to facilitate mRNA biosynthesis, transcriptional regulation, and signal transduction. Research studies have shown that cytoplasmic hnRNP K expression is increased in oral squamous cell carcinoma and pancreatic cancer, and may be a potential prognostic factor.


Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of hnRNP K in K562, rat Brain, C6, 3T3, HeLa lysates using hnRNP K antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using hnRNP K antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.