

Product Name: HPDL Rabbit Monoclonal Antibody**Catalog #: AMRe02116**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:20-1:50
Molecular Weight	Calculated MW: 39 kDa; Observed MW: 39 kDa

Antigen Information

Gene Name	HPDL
Alternative Names	4-hydroxyphenylpyruvate dioxygenase like; GLOXD1; 4-HPPD-L
Gene ID	84842
SwissProt ID	Q96IR7
Immunogen	A synthetic peptide of human HPDL


Background

May have dioxygenase activity.

Research Area

Signal Transduction

Image Data



Western blot analysis of HPDL in HeLa, Caco2, A549 lysates using HPDL antibody.