

Product Name: INPP4B Rabbit Monoclonal Antibody
Catalog #: AMRe02161



Summary

Production Name	INPP4B Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	INPP4B 4-bisphosphate 4-phosphatase; Inositol polyphosphate 4 phosphatase type II 105kDa;
Alternative Names	Inositol polyphosphate 4-phosphatase type II; Inpp4b; MGC132014; Type II inositol 3; 4 bisphosphate 4 phosphatase; Type II inositol-3
Gene ID	8821
SwissProt ID	O15327.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 105 kDa; Observed MW: 105 kDa

Product Name: INPP4B Rabbit Monoclonal Antibody
Catalog #: AMRe02161




Background

Catalyzes the hydrolysis of the 4-position phosphate of phosphatidylinositol 3,4-bisphosphate, inositol 1,3,4-trisphosphate and inositol 1,4-bisphosphate.


Research Area

Signal Transduction

Image Data



Western blot analysis of INPP4B in Hacat lysates using INPP4B antibody.



Western blot analysis of INPP4B in rat Brain, C6 lysates using INPP4B antibody.

Note

For research use only.