
Product Name: IRF8 Rabbit Monoclonal Antibody**Catalog #: AMRe02172**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
Molecular Weight	Calculated MW: 48 kDa; Observed MW: 48 kDa

Antigen Information

Gene Name	IRF8
Alternative Names	ICSBP; IRF-8; ICSBP1; IMD32A; IMD32B; H-ICSBP
Gene ID	3394
SwissProt ID	Q02556
Immunogen	Recombinant protein of human IRF8

Background


Plays a role as a transcriptional activator or repressor (PubMed:25122610). Specifically binds to the upstream regulatory region of type I IFN and IFN-inducible MHC class I genes (the interferon consensus sequence (ICS)). Plays a negative regulatory role in

cells of the immune system. Involved in CD8+ dendritic cell differentiation by forming a complex with the BATF-JUNB heterodimer in immune cells, leading to recognition of AICE sequence (5'-TGAnTCA/GAAA-3'), an immune-specific regulatory element, followed by cooperative binding of BATF and IRF8 and activation of genes .

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of IRF8 in THP-1 lysates using IRF8 antibody.