

---

**Product Name: ISG15 Rabbit Monoclonal antibody****Catalog #: AMRe02173**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal Antibody
<b>Form</b>	Liquid
<b>Concentration</b>	0.5mg/ml. The concentration of this product may be batch-dependent.
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
<b>Purification</b>	Affinity Purified

**Application**

<b>Dilution Ratio</b>	WB 1:500-1:1000
<b>Molecular Weight</b>	Calculated MW: 18 kDa; Observed MW: 18 kDa

**Antigen Information**

<b>Gene Name</b>	ISG15
<b>Alternative Names</b>	IP17; Ubiquitin cross-reactive protein; hUCRP; ISG15; G1P2
<b>Gene ID</b>	9636
<b>SwissProt ID</b>	P05161
<b>Immunogen</b>	A synthetic peptide of human ISG15

**Background**


May serve as a trans-acting binding factor directing the association of ligated target proteins to intermediate filaments. May also be involved in autocrine, paracrine and endocrine mechanisms, as in cell-to-cell signaling, possibly partly by inducing IFN-

gamma secretion by monocytes and macrophages. Seems to display antiviral activity during viral infections.

## Research Area

Cell Biology

## Image Data



Western blot analysis of ISG15 in C6, 3T3, HeLa lysates using ISG15 antibody.