

Product Name: KDM4B Rabbit Monoclonal Antibody
Catalog #: AMRe02191



Summary

Production Name	KDM4B Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	KDM4B
Alternative Names	KDM4B; JHDM3B; JMJD2B; KIAA0876; Lysine-specific demethylase 4B; JmjC domain-containing histone demethylation protein 3B; Jumonji domain-containing protein 2B
Gene ID	23030
SwissProt ID	O94953.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IF: 1:50-1:200
Molecular Weight	Calculated MW: 122 kDa; Observed MW: 150 kDa

Product Name: KDM4B Rabbit Monoclonal Antibody
Catalog #: AMRe02191




Background

Histone demethylase that specifically demethylates 'Lys-9' of histone H3, thereby playing a role in histone code. Does not demethylate histone H3 'Lys-4', H3 'Lys-27', H3 'Lys-36' nor H4 'Lys-20'. Only able to demethylate trimethylated H3 'Lys-9', with a weaker activity than KDM4A, KDM4C and KDM4D. Demethylation of Lys residue generates formaldehyde and succinate.


Research Area

Epigenetics and Nuclear Signaling

Image Data




Western blot analysis of KDM4B in 293, COS7 lysates using KDM4B antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using KDM4B/JMJD2B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Product Name: KDM4B Rabbit Monoclonal Antibody
Catalog #: AMRe02191



Immunocytochemistry analysis of KDM4B/JMJD2B(green) in HeLa using KDM4B/JMJD2B antibody, and DAPI(blue)

Note

For research use only.