
Product Name: LGALS3BP Rabbit Monoclonal Antibody**Catalog #: AMRe02213**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	1.18mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 65 kDa; Observed MW: 100 kDa

Antigen Information

Gene Name	LGALS3BP
Alternative Names	90K; M2BP; gp90; CyCAP; BTBD17B; MAC-2-BP; TANGO10B
Gene ID	3959
SwissProt ID	Q08380
Immunogen	Recombinant protein of human LGALS3BP


Background

Promotes integrin-mediated cell adhesion. May stimulate host defense against viruses and tumor cells.

Research Area

Cancer

Image Data



Western blot analysis of LGALS3BP in HT-1080, HepG2, A549 lysates using LGALS3BP antibody.