
Product Name: LRG1 Rabbit Monoclonal Antibody**Catalog #: AMRe02220**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
Molecular Weight	Calculated MW: 38 kDa; Observed MW: 48 kDa

Antigen Information

Gene Name	LRG1
Alternative Names	LRG1; LRG; Leucine-rich alpha-2-glycoprotein; LRG
Gene ID	116844
SwissProt ID	P02750
Immunogen	A synthetic peptide of human LRG1

Background

Transforming growth factor beta receptor binding.

Research Area

Immunology

Image Data

Western blot analysis of LRG1 in Human serum lysates using LRG1 antibody.

