

**Product Name: MCM2 Rabbit Monoclonal antibody**  
**Catalog #: AMRe02234**



## Summary

|                        |                                        |
|------------------------|----------------------------------------|
| <b>Production Name</b> | MCM2 Rabbit Monoclonal antibody        |
| <b>Description</b>     | Recombinant rabbit monoclonal antibody |
| <b>Host</b>            | Rabbit                                 |
| <b>Application</b>     | WB, ICC/IF                             |
| <b>Reactivity</b>      | Human, Mouse, Rat                      |

## Performance

|                     |                                                                                                      |
|---------------------|------------------------------------------------------------------------------------------------------|
| <b>Conjugation</b>  | Unconjugated                                                                                         |
| <b>Modification</b> | Unmodified                                                                                           |
| <b>Isotype</b>      | IgG                                                                                                  |
| <b>Clonality</b>    | Monoclonal Antibody                                                                                  |
| <b>Form</b>         | Liquid                                                                                               |
| <b>Storage</b>      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.             |
| <b>Buffer</b>       | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein |
| <b>Purification</b> | Affinity Purified                                                                                    |

## Immunogen

|                          |                                                                                                                                               |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | MCM2                                                                                                                                          |
| <b>Alternative Names</b> | MCM2; BM28; CCNL1; CDCL1; KIAA0030; DNA replication licensing factor MCM2; Minichromosome maintenance protein 2 homolog; Nuclear protein BM28 |
| <b>Gene ID</b>           | 4171                                                                                                                                          |
| <b>SwissProt ID</b>      | P49736.                                                                                                                                       |

## Application

|                         |                                              |
|-------------------------|----------------------------------------------|
| <b>Dilution Ratio</b>   | WB 1:500-1:1000, ICC/IF 1:50-1:200           |
| <b>Molecular Weight</b> | Calculated MW: 102 kDa; Observed MW: 125 kDa |

**Product Name: MCM2 Rabbit Monoclonal antibody**  
**Catalog #: AMRe02234**




## Background

The protein encoded by this gene is one of the highly conserved mini-chromosome maintenance proteins (MCM) that are involved in the initiation of eukaryotic genome replication. The hexameric protein complex formed by MCM proteins is a key component of the pre-replication complex (pre\_RC) and may be involved in the formation of replication forks and in the recruitment of other DNA replication related proteins. This protein forms a complex with MCM4, 6, and 7, and has been shown to regulate the helicase activity of the complex. This protein is phosphorylated, and thus regulated by, protein kinases CDC2 and CDC7. Multiple alternatively spliced transcript variants have been found, but the full-length nature of some variants has not been defined. [provided by RefSeq, Oct 2012]


## Research Area

Epigenetics and Nuclear Signaling

## Image Data



Western blot analysis of MCM2 in K562, C6, 3T3, HeLa lysates using MCM2 antibody.



Immunocytochemistry analysis of MCM2(green) in HeLa using MCM2 antibody, and DAPI(blue)

## Note

For research use only.