

---

**Product Name: MTCO2 Rabbit Monoclonal Antibody****Catalog #: AMRe02280**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC,ICC/IF,IP
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	0.61mg/ml. The concentration of this product may be batch-dependent.
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
<b>Molecular Weight</b>	Calculated MW: 26 kDa; Observed MW: 21 kDa

**Antigen Information**

<b>Gene Name</b>	MT-CO2
<b>Alternative Names</b>	MT-CO2; COX2; CO2; Cytochrome c oxidase II; MTCO2; COII; COXII
<b>Gene ID</b>	4513
<b>SwissProt ID</b>	P00403
<b>Immunogen</b>	A synthetic peptide of human MTCO2

**Background**


Cytochrome c oxidase is the component of the respiratory chain that catalyzes the reduction of oxygen to water. Subunits 1-3 form the functional core of the enzyme complex. Subunit 2 transfers the electrons from cytochrome c via its binuclear copper A

center to the bimetallic center of the catalytic subunit 1.


## Research Area

Signal Transduction


## Image Data



Western blot analysis of MTCO2 in K562, HeLa lysates using MTCO2 antibody.



Immunocytochemistry analysis of MTCO2(green) in HeLa using MTCO2 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using MTCO2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.