

**Product Name: Muscarinic Acetylcholine Receptor M2
Rabbit Monoclonal Antibody
Catalog #: AMRe02290**

Summary

Production Name	Muscarinic Acetylcholine Receptor M2 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P,IP
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	CHRM2
Alternative Names	CHRM2; HM2; AChR; Acm2
Gene ID	1129
SwissProt ID	P08172.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IP: 1:20
Molecular Weight	Calculated MW: 52 kDa; Observed MW: 52 kDa

Background

**Product Name: Muscarinic Acetylcholine Receptor M2
Rabbit Monoclonal Antibody
Catalog #: AMRe02290**




The muscarinic acetylcholine receptor mediates various cellular responses, including inhibition of adenylate cyclase, breakdown of phosphoinositides and modulation of potassium channels through the action of G proteins. Primary transducing effect is adenylate cyclase inhibition.

Research Area

Neuroscience

Image Data



Western blot analysis of Muscarinic Acetylcholine Receptor 2/CM2 in K562, C6, HeLa lysates using Muscarinic Acetylcholine Receptor M2 antibody.

Note

For research use only.