

Product Name: MVP Rabbit Monoclonal Antibody
Catalog #: AMRe02292



Summary

Production Name	MVP Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	MVP
Alternative Names	LRP; MVP; VAULT1
Gene ID	9961
SwissProt ID	Q14764.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IF: 1:50-1:200
Molecular Weight	Calculated MW: 99 kDa; Observed MW: 110 kDa

Background

Required for normal vault structure. Vaults are multi-subunit structures that may act as scaffolds for proteins involved in

Product Name: MVP Rabbit Monoclonal Antibody
Catalog #: AMRe02292




signal transduction. Vaults may also play a role in nucleo-cytoplasmic transport.

Research Area

Signal Transduction

Image Data



Western blot analysis of MVP in U87-MG lysates using MVP antibody.

Note

For research use only.