

---

**Product Name: NFIB Rabbit Monoclonal antibody****Catalog #: AMRe02340**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IP
<b>Reactivity</b>	Human,Rat
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal Antibody
<b>Form</b>	Liquid
<b>Concentration</b>	1.47mg/ml. The concentration of this product may be batch-dependent.
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
<b>Purification</b>	Affinity Purified

**Application**

<b>Dilution Ratio</b>	WB 1:500-1:1000,IP 1:20-1:50
<b>Molecular Weight</b>	Calculated MW: 47 kDa; Observed MW: 47,55 kDa

**Antigen Information**

<b>Gene Name</b>	NFIB
<b>Alternative Names</b>	CTF; HMGIC/NFIB; NF-I/B; NF1-B; NF1B; NF1B2; NFI-B; NFI-RED; Nfib; NFIB2; NFIB3
<b>Gene ID</b>	4781
<b>SwissProt ID</b>	O00712
<b>Immunogen</b>	A synthetic peptide of human NFIB


**Background**

Recognizes and binds the palindromic sequence 5'-TTGGCNNNNNGCCAA-3' present in viral and cellular promoters and in the origin of replication of adenovirus type 2. These proteins are individually capable of activating transcription and replication.

## Research Area

Neuroscience

## Image Data



Western blot analysis of NFIB in C6, HeLa lysates using NFIB antibody.