

Product Name: Nova1 Rabbit Monoclonal antibody**Catalog #: AMRe02352**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Concentration	1.47mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100
Molecular Weight	Calculated MW: 52 kDa; Observed MW: 52 kDa

Antigen Information

Gene Name	NOVA1
Alternative Names	Nova-1
Gene ID	4857
SwissProt ID	P51513
Immunogen	A synthetic peptide of human Nova1


Background

May regulate RNA splicing or metabolism in a specific subset of developing neurons.

Research Area

Neuroscience

Image Data



Western blot analysis of Nova1 in rat lung lysates using Nova1 antibody.