

Product Name: NTH1 Rabbit Monoclonal antibody
Catalog #: AMRe02358



Summary

Production Name	NTH1 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	NTHL1
Alternative Names	FAP3; NTH1; OCTS3; hNTH1
Gene ID	4913
SwissProt ID	P78549.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200
Molecular Weight	Calculated MW: 34 kDa; Observed MW: 34 kDa

Background

Product Name: NTH1 Rabbit Monoclonal antibody
Catalog #: AMRe02358




Bifunctional DNA N-glycosylase with associated apurinic/aprimidinic (AP) lyase function that catalyzes the first step in base excision repair (BER), the primary repair pathway for the repair of oxidative DNA damage. The DNA N-glycosylase activity releases the damaged DNA base from DNA by cleaving the N-glycosidic bond, leaving an AP site. The AP-lyase activity cleaves the phosphodiester bond 3' to the AP site by a beta-elimination. Primarily recognizes and repairs oxidative base damage of pyrimidines. Has also 8-oxo-7,8-dihydroguanine (8-oxoG) DNA glycosylase activity. Acts preferentially on DNA damage opposite guanine residues in DNA. Is able to process lesions in nucleosomes without requiring or inducing nucleosome disruption.


Research Area

Epigenetics and Nuclear Signaling

Image Data



Immunocytochemistry analysis of NTH1 (green) in HeLa using NTH1 antibody, and DAPI (blue).



Western blot analysis of NTH1 in K562, HeLa lysates using NTH1 antibody.

Note

For research use only.