

---

**Product Name: NXF1 Rabbit Monoclonal Antibody****Catalog #: AMRe02370**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC,ICC/IF,IP
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	0.54mg/ml. The concentration of this product may be batch-dependent.
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
<b>Molecular Weight</b>	Calculated MW: 70 kDa; Observed MW: 70 kDa

**Antigen Information**

<b>Gene Name</b>	NXF1
<b>Alternative Names</b>	DmNXF1; MEX67; Mex67p; Mvb1; NXF1; Sbr; TAP
<b>Gene ID</b>	10482
<b>SwissProt ID</b>	Q9UBU9
<b>Immunogen</b>	A synthetic peptide of human NXF1

**Background**


Involved in the nuclear export of mRNA species bearing retroviral constitutive transport elements (CTE) and in the export of mRNA from the nucleus to the cytoplasm. The NXF1-NXT1 heterodimer is involved in the export of HSP70 mRNA in

conjunction with THOC4 and THOC5.


## Research Area

Tags & Cell Markers


## Image Data



Immunocytochemistry analysis of NXF1 (green) in Jurkat using NXF1 antibody, and DAPI (blue).



Western blot analysis of NXF1 in K562 lysates using NXF1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer tissue using NXF1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. A 40µm scale bar is present in the bottom right corner.