

Product Name: P Glycoprotein Rabbit Monoclonal antibody
Catalog #: AMRe02378

Summary

Production Name	P Glycoprotein Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	ABCB1
Alternative Names	p-pg; ABCB1; MDR1; PGY1; Multidrug resistance protein 1; ATP-binding cassette sub-family B member 1; P-glycoprotein 1; CD antigen CD243
Gene ID	5243
SwissProt ID	P08183.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 141 kDa; Observed MW: 130-180 kDa

Product Name: P Glycoprotein Rabbit Monoclonal antibody
Catalog #: AMRe02378


Background

A chronic, relapsing inflammation of the gastrointestinal tract with a complex etiology. It is subdivided into Crohn disease and ulcerative colitis phenotypes. Crohn disease may affect any part of the gastrointestinal tract from the mouth to the anus, but most frequently it involves the terminal ileum and colon.

Research Area

Signal Transduction

Image Data



Western blot analysis of P Glycoprotein in HeLa lysates using P Glycoprotein antibody.

Note

For research use only.