

Product Name: PDZK1 Rabbit Monoclonal Antibody**Catalog #: AMRe02418**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:20-1:50
Molecular Weight	Calculated MW: 57 kDa; Observed MW: 70 kDa

Antigen Information

Gene Name	PDZK1
Alternative Names	CAP70; CLAMP; PDZD1; NHERF3; NHERF-3
Gene ID	5174
SwissProt ID	Q5T2W1
Immunogen	A synthetic peptide of human PDZK1


Background

A scaffold protein that connects plasma membrane proteins and regulatory components, regulating their surface expression in epithelial cells apical domains.

Research Area

Signal Transduction

Image Data



Western blot analysis of PDZK1 in T4D7, Caco2 lysates using PDZK1 antibody.