

Product Name: PHLDA1 Rabbit Monoclonal Antibody
Catalog #: AMRe02435



Summary

Production Name	PHLDA1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	PHLDA1
Alternative Names	PQ-rich protein; PQR protein; Proline- and histidine-rich protein; T-cell death-associated gene 51 protein
Gene ID	22822
SwissProt ID	Q8WV24.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IF: 1:50-1:200
Molecular Weight	Calculated MW: 45 kDa; Observed MW: 45 kDa

Background

Product Name: PHLDA1 Rabbit Monoclonal Antibody
Catalog #: AMRe02435




Seems to be involved in regulation of apoptosis.

Research Area

Cell Biology

Image Data



Western blot analysis of PHLDA1 in A375 lysates using PHLDA1 antibody.

Note

For research use only.