

Product Name: PIN4 Rabbit Monoclonal Antibody
Catalog #: AMRe02443



Summary

Production Name	PIN4 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Rat,Hamster

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	PIN4
Alternative Names	EPVH; PAR14; PAR17
Gene ID	5303
SwissProt ID	Q9Y237.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 14 kDa; Observed MW: 14 kDa

Background

Isoform 1 is involved as a ribosomal RNA processing factor in ribosome biogenesis. Binds to tightly bent AT-rich stretches

Product Name: PIN4 Rabbit Monoclonal Antibody
Catalog #: AMRe02443




of double-stranded DNA.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of PIN4 in CHO-K1, C6, Jurkat lysates using PIN4 antibody.

Note

For research use only.