

---

**Product Name: PON1 Rabbit Monoclonal Antibody****Catalog #: AMRe02458**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IP
<b>Reactivity</b>	Human,Mouse,Rat
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	0.67mg/ml. The concentration of this product may be batch-dependent.
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	WB 1:500-1:1000,IP 1:20-1:50
<b>Molecular Weight</b>	Calculated MW: 40 kDa; Observed MW: 40 kDa

**Antigen Information**

<b>Gene Name</b>	PON1
<b>Alternative Names</b>	A esterase 1; Arylesterase 1; ESA; MVCD5; Paraoxonase 1; PON1
<b>Gene ID</b>	5444
<b>SwissProt ID</b>	P27169
<b>Immunogen</b>	A synthetic peptide of human PON1

**Background**


Hydrolyzes the toxic metabolites of a variety of organophosphorus insecticides. Capable of hydrolyzing a broad spectrum of organophosphate substrates and lactones, and a number of aromatic carboxylic acid esters. Mediates an enzymatic protection

of low density lipoproteins against oxidative modification and the consequent series of events leading to atheroma formation.


## Research Area

Neuroscience

## Image Data



Western blot analysis of PON1 in mouse liver lysates using PON1 antibody.



Western blot analysis of PON1 in K562 lysates using PON1 antibody.