

Product Name: Prealbumin Rabbit Monoclonal Antibody
Catalog #: AMRe02470

Summary

Production Name	Prealbumin Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	TTR
Alternative Names	Prealbumin; TBPA; Amyloidosis I; Transthyretin
Gene ID	7276
SwissProt ID	P02766.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 16 kDa; Observed MW: 16 kDa

Background

Thyroid hormone-binding protein. Probably transports thyroxine from the bloodstream to the brain.


Product Name: Prealbumin Rabbit Monoclonal Antibody
Catalog #: AMRe02470



Research Area

Neuroscience

Image Data



Western blot analysis of Prealbumin in rat Brain lysates using Prealbumin antibody.

Note

For research use only.