

Product Name: Prealbumin Rabbit Monoclonal Antibody
Catalog #: AMRe02471

Summary

Production Name	Prealbumin Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	TTR
Alternative Names	Prealbumin; TBPA; Amyloidosis I; Transthyretin
Gene ID	7276
SwissProt ID	P02766.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IP: 1:20
Molecular Weight	Calculated MW: 16 kDa; Observed MW: 16 kDa

Background


Thyroid hormone-binding protein. Probably transports thyroxine from the bloodstream to the brain.

Product Name: Prealbumin Rabbit Monoclonal Antibody
Catalog #: AMRe02471

Research Area

Neuroscience

Image Data



Western blot analysis of Prealbumin in Human tonsil, Human placenta lysates using Prealbumin antibody.

Note

For research use only.