

Product Name: PKA C gamma Rabbit Monoclonal Antibody
Catalog #: AMRe02473

Summary

Production Name	PKA C gamma Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	PRKACG
Alternative Names	PRKACG; cAMP-dependent protein kinase catalytic subunit gamma; PKA C-gamma
Gene ID	5568
SwissProt ID	P22612.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 40 kDa; Observed MW: 40 kDa

Product Name: PKA C gamma Rabbit Monoclonal Antibody
Catalog #: AMRe02473




Background

Phosphorylates a large number of substrates in the cytoplasm and the nucleus.


Research Area

Signal Transduction

Image Data



Western blot analysis of PRKACG in mouse testis, mouse lung lysates using PKA C gamma antibody.



Western blot analysis of PRKACG in MCF-7, C6, HeLa lysates using PRKACG antibody.

Note

For research use only.