

Product Name: PRMT3 Rabbit Monoclonal Antibody**Catalog #: AMRe02474**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.67mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
Molecular Weight	Calculated MW: 60 kDa; Observed MW: 60 kDa

Antigen Information

Gene Name	PRMT3
Alternative Names	HRMT1L3
Gene ID	10196
SwissProt ID	O60678
Immunogen	A synthetic peptide of human PRMT3


Background

Methylates (mono and asymmetric dimethylation) the guanidino nitrogens of arginyl residues in some proteins.


Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of PRMT3 in rat Brain, C6, 3T3, Hela lysates using PRMT3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using PRMT3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.