

Product Name: PTN Rabbit Monoclonal antibody
Catalog #: AMRe02495



Summary

Production Name	PTN Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	PTN
Alternative Names	PTN; HBNF1; NEGF1; Pleiotrophin; PTN; Heparin-binding brain mitogen; HBBM; Heparin-binding growth factor 8; HBGF-8; Heparin-binding growth-associated molecule; HB-GAM; Heparin-binding neurite outgrowth-promoting factor 1; HBNF-1; Osteoblast-specific factor 1; OSF-1
Gene ID	5764
SwissProt ID	P21246.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 19 kDa; Observed MW: 19 kDa

Product Name: PTN Rabbit Monoclonal antibody
Catalog #: AMRe02495




Background

Secreted growth factor that induces neurite outgrowth and which is mitogenic for fibroblasts, epithelial, and endothelial cells (PubMed:1768439, PubMed:1733956). Binds anaplastic lymphoma kinase (ALK) which induces MAPK pathway activation, an important step in the anti-apoptotic signaling of PTN and regulation of cell proliferation (PubMed:11278720). Binds to cell-surface target proteins via their chondroitin sulfate groups (PubMed:26896299). Down-regulates PTPRZ1 activity (PubMed:16814777).


Research Area

Neuroscience

Image Data



Western blot analysis of PTN in U87-MG lysates using PTN antibody.



Western blot analysis of PTN in C6, rat Brain lysates using PTN antibody.

Note

For research use only.

Product Name: PTN Rabbit Monoclonal antibody
Catalog #: AMRe02495

