
Product Name: PWP1 Rabbit Monoclonal antibody**Catalog #: AMRe02501**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 56 kDa; Observed MW: 75 kDa

Antigen Information

Gene Name	PWP1
Alternative Names	IEF-SSP-9502
Gene ID	11137
SwissProt ID	Q13610
Immunogen	A synthetic peptide of human PWP1


Background

Chromatin-associated factor that regulates transcription (PubMed:29065309). Regulates Pol I-mediated rRNA biogenesis and, probably, Pol III-mediated transcription (PubMed:29065309). Regulates the epigenetic status of rDNA (PubMed:29065309).

Research Area

Cell Biology

Image Data



Western blot analysis of PWP1 in K562 lysates using PWP1 antibody.