

Product Name: Rab22A Rabbit Monoclonal antibody
Catalog #: AMRe02505



Summary

Production Name	Rab22A Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	RAB22A
Alternative Names	Ras-related protein Rab-22A
Gene ID	57403
SwissProt ID	Q9UL26.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 22 kDa; Observed MW: 22 kDa

Background

Plays a role in endocytosis and intracellular protein transport. Mediates trafficking of TF from early endosomes to recycling

Product Name: Rab22A Rabbit Monoclonal antibody
Catalog #: AMRe02505




endosomes (PubMed:16537905). Required for NGF-mediated endocytosis of NTRK1, and subsequent neurite outgrowth (PubMed:21849477). Binds GTP and GDP and has low GTPase activity. Alternates between a GTP-bound active form and a GDP-bound inactive form (PubMed:16537905).

Research Area

Signal Transduction

Image Data



Western blot analysis of RAB22A in rat Brain, C6, 3T3, HeLa lysates using Rab22A antibody.

Note

For research use only.