

Product Name: Rab9 Rabbit Monoclonal Antibody
Catalog #: AMRe02510



Summary

Production Name	Rab9 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	RAB9A
Alternative Names	RAB9A; DmRab9; Ras related protein Rab 9A
Gene ID	9367
SwissProt ID	P51151.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100
Molecular Weight	Calculated MW: 23 kDa; Observed MW: 23 kDa

Background

Product Name: Rab9 Rabbit Monoclonal Antibody
Catalog #: AMRe02510




Involved in the transport of proteins between the endosomes and the trans Golgi network. Interacts (preferentially in its GTP-bound form) with GCC2 (via its GRIP domain).


Research Area

Tags & Cell Markers

Image Data



Western blot analysis of Rab9 in 3T3, HeLa lysates using Rab9 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Rab9 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.