

Product Name: REEP2 Rabbit Monoclonal Antibody
Catalog #: AMRe02535



Summary

Production Name	REEP2 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	REEP2
Alternative Names	Receptor expression-enhancing protein 2
Gene ID	51308
SwissProt ID	Q9BRK0.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 28 kDa; Observed MW: 28 kDa

Background

Required for endoplasmic reticulum (ER) network formation, shaping and remodeling.


Product Name: REEP2 Rabbit Monoclonal Antibody
Catalog #: AMRe02535



Research Area

Neuroscience

Image Data



Western blot analysis of REEP2 in Raoms, Jurkat lysates using REEP2 antibody.

Note

For research use only.