

Product Name: Sall4 Rabbit Monoclonal Antibody**Catalog #: AMRe02570**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.3mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100
Molecular Weight	Calculated MW: 112 kDa; Observed MW: 112 kDa

Antigen Information

Gene Name	SALL4
Alternative Names	Zinc finger protein 797; Zinc finger protein SALL4
Gene ID	57167
SwissProt ID	Q9UJQ4
Immunogen	Recombinant protein of human Sall4


Background

Transcription factor with a key role in the maintenance and self-renewal of embryonic and hematopoietic stem cells.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of Sall4 in Daudi lysates using Sall4 antibody.