

Product Name: Semaphorin 4D Rabbit Monoclonal antibody
Catalog #: AMRe02585

Summary

Production Name	Semaphorin 4D Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	SEMA4D
Alternative Names	CD100; COLL4; SEMAJ; coll-4; C9orf164; M-sema-G
Gene ID	10507
SwissProt ID	Q92854.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200
Molecular Weight	Calculated MW: 96 kDa; Observed MW: 120 kDa

Product Name: Semaphorin 4D Rabbit Monoclonal antibody
Catalog #: AMRe02585




Background

Cell surface receptor for PLXN1B and PLXNB2 that plays an important role in cell-cell signaling. Promotes reorganization of the actin cytoskeleton and plays a role in axonal growth cone guidance in the developing central nervous system. Regulates dendrite and axon branching and morphogenesis. Promotes the migration of cerebellar granule cells and of endothelial cells. Plays a role in the immune system; induces B-cells to aggregate and improves their viability (in vitro). Promotes signaling via SRC and PTK2B/PYK2, which then mediates activation of phosphatidylinositol 3-kinase and of the AKT1 signaling cascade. Interaction with PLXNB1 mediates activation of RHOA.


Research Area

Immunology

Image Data



Immunocytochemistry analysis of Semaphorin 4D (green) in Jurkat using Semaphorin 4D antibody, and DAPI (blue).



Western blot analysis of Semaphorin 4D/CD100 in Jurkat lysates using Semaphorin 4D/CD100 antibody.

Note

For research use only.