

Product Name: SLCO5A1 Rabbit Monoclonal antibody
Catalog #: AMRe02611



Summary

Production Name	SLCO5A1 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	SLCO5A1
Alternative Names	OATPJ; OATP-J; OATP5A1; OATPRP4; OATP-RP4; SLC21A15
Gene ID	81796
SwissProt ID	Q9H2Y9.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100
Molecular Weight	Calculated MW: 92 kDa; Observed MW: 92 kDa

Background

Product Name: SLCO5A1 Rabbit Monoclonal antibody
Catalog #: AMRe02611




Solute carrier organic anion transporter family member 5A1


Research Area

Signal Transduction

Image Data



Western blot analysis of SLCO5A1 in T47D lysates using SLCO5A1 antibody.



Western blot analysis of SLCO5A1 in K562, HeLa lysates using SLCO5A1 antibody.

Note

For research use only.