

Product Name: SNX4 Rabbit Monoclonal antibody
Catalog #: AMRe02624



Summary

Production Name	SNX4 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	SNX4
Alternative Names	ATG24B
Gene ID	8723
SwissProt ID	O95219.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 52 kDa; Observed MW: 52 kDa

Background

Involved in the regulation of endocytosis and in several stages of intracellular trafficking (PubMed:12668730). Plays a role

Product Name: SNX4 Rabbit Monoclonal antibody
Catalog #: AMRe02624




in recycling endocytosed transferrin receptor and prevent its degradation.

Research Area

Signal Transduction

Image Data



Western blot analysis of SNX4 in C6, 3T3, Hela lysates using SNX4 antibody.

Note

For research use only.