

Product Name: SPOP Rabbit Monoclonal antibody
Catalog #: AMRe02633



Summary

Production Name	SPOP Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	SPOP
Alternative Names	Speckle-type POZ protein; HIB homolog 1; Roadkill homolog 1
Gene ID	8405
SwissProt ID	O43791.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 42 kDa; Observed MW: 42 kDa

Background

Inhibits IPF1/PDX1 transactivation of established target promoters, such as insulin, may be by recruiting a repressor

Product Name: SPOP Rabbit Monoclonal antibody
Catalog #: AMRe02633




complex (By similarity). In complex with CUL3, involved in ubiquitination of BMI1, H2AFY and DAXX, and probably also in ubiquitination and proteasomal degradation of Gli2 or Gli3.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of SPOP in MCF-7 lysates using SPOP antibody.

Note

For research use only.