

**Product Name: STIP1 Rabbit Monoclonal antibody**  
**Catalog #: AMRe02642**



## Summary

|                        |                                        |
|------------------------|----------------------------------------|
| <b>Production Name</b> | STIP1 Rabbit Monoclonal antibody       |
| <b>Description</b>     | Recombinant Rabbit Monoclonal antibody |
| <b>Host</b>            | Rabbit                                 |
| <b>Application</b>     | WB,IHC-P,IP                            |
| <b>Reactivity</b>      | Human,Mouse,Rat                        |

## Performance

|                     |                                                                                                      |
|---------------------|------------------------------------------------------------------------------------------------------|
| <b>Conjugation</b>  | Unconjugated                                                                                         |
| <b>Modification</b> | Unmodified                                                                                           |
| <b>Isotype</b>      | IgG                                                                                                  |
| <b>Clonality</b>    | Monoclonal Antibody                                                                                  |
| <b>Form</b>         | Liquid                                                                                               |
| <b>Storage</b>      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.             |
| <b>Buffer</b>       | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein |
| <b>Purification</b> | Affinity Purified                                                                                    |

## Immunogen

|                          |                                                               |
|--------------------------|---------------------------------------------------------------|
| <b>Gene Name</b>         | STIP1                                                         |
| <b>Alternative Names</b> | Hop; P60; STI1; STI1L; STIP1; Stress induced phosphoprotein 1 |
| <b>Gene ID</b>           | 10963                                                         |
| <b>SwissProt ID</b>      | P31948.                                                       |

## Application

|                         |                                            |
|-------------------------|--------------------------------------------|
| <b>Dilution Ratio</b>   | WB: 1:500-1:1000 IHC: 1:50-1:100 IP: 1:20  |
| <b>Molecular Weight</b> | Calculated MW: 63 kDa; Observed MW: 63 kDa |

## Background

**Product Name: STIP1 Rabbit Monoclonal antibody**  
**Catalog #: AMRe02642**




Mediates the association of the molecular chaperones HSC70 and HSP90 (HSPCA and HSPCB).


## Research Area

Signal Transduction


## Image Data



Western blot analysis of STIP1 in mouse brain, mouse heart lysates using STIP1 antibody.



Western blot analysis of STIP1 in K562, rat Brain, C6 lysates using STIP1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using STIP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

**Product Name: STIP1 Rabbit Monoclonal antibody**  
**Catalog #: AMRe02642**

---



**Note**

For research use only.