

Product Name: Syndecan 3 Rabbit Monoclonal Antibody
Catalog #: AMRe02656

Summary

Production Name	Syndecan 3 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	SDC3
Alternative Names	SDCN; SYND3
Gene ID	9672
SwissProt ID	O75056.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 45 kDa; Observed MW: 60 kDa

Background

Cell surface proteoglycan that may bear heparan sulfate . May have a role in the organization of cell shape by affecting the

Product Name: Syndecan 3 Rabbit Monoclonal Antibody
Catalog #: AMRe02656




actin cytoskeleton, possibly by transferring signals from the cell surface in a sugar-dependent mechanism.

Research Area

Signal Transduction

Image Data



Western blot analysis of Syndecan 3 in K562, C6, 3T3 lysates using Syndecan 3 antibody.

Note

For research use only.