

Product Name: Syntaxin 16 Rabbit Monoclonal antibody**Catalog #: AMRe02657**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Concentration	0.28mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 37 kDa; Observed MW: 37 kDa

Antigen Information

Gene Name	STX16
Alternative Names	syntaxin 16; SYN16
Gene ID	8675
SwissProt ID	O14662
Immunogen	A synthetic peptide of human Syntaxin 16


Background

SNARE involved in vesicular transport from the late endosomes to the trans-Golgi network.


Research Area

Tags & Cell Markers

Image Data



Western blot analysis of Syntaxin 16 in mouse lung lysates using Syntaxin 16 antibody.



Western blot analysis of Syntaxin 16 in K562, rat Brain lysates using Syntaxin 16 antibody.