

Product Name: TAGAP Rabbit Monoclonal antibody**Catalog #: AMRe02663**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Concentration	1.68mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:20-1:50
Molecular Weight	Calculated MW: 81 kDa; Observed MW: 81 kDa

Antigen Information

Gene Name	TAGAP
Alternative Names	FKSG15; IDDM21; TAGAP1; ARHGAP47
Gene ID	117289
SwissProt ID	Q8N103
Immunogen	Recombinant protein of human TAGAP


Background

May function as a GTPase-activating protein and may play important roles during T-cell activation.

Research Area

Signal Transduction

Image Data



Western blot analysis of TAGAP in rat Brain, C6, CHO-K1, Hela lysates using TAGAP antibody.