

**Product Name: T Box Protein 3 Rabbit Monoclonal Antibody**  
**Catalog #: AMRe02671**

---

## Summary

|                        |                                            |
|------------------------|--------------------------------------------|
| <b>Production Name</b> | T Box Protein 3 Rabbit Monoclonal Antibody |
| <b>Description</b>     | Rabbit Monoclonal antibody                 |
| <b>Host</b>            | Rabbit                                     |
| <b>Application</b>     | WB,IP                                      |
| <b>Reactivity</b>      | Mouse                                      |

## Performance

|                     |                                                                                                      |
|---------------------|------------------------------------------------------------------------------------------------------|
| <b>Conjugation</b>  | Unconjugated                                                                                         |
| <b>Modification</b> | Unmodified                                                                                           |
| <b>Isotype</b>      | IgG                                                                                                  |
| <b>Clonality</b>    | Monoclonal                                                                                           |
| <b>Form</b>         | Liquid                                                                                               |
| <b>Storage</b>      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.             |
| <b>Buffer</b>       | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein |
| <b>Purification</b> | Affinity Purification                                                                                |

## Immunogen

|                          |                                                    |
|--------------------------|----------------------------------------------------|
| <b>Gene Name</b>         | TBX3                                               |
| <b>Alternative Names</b> | T box 3; T-box protein 3; Tbx3; TBX3 ISO; UMS; XHL |
| <b>Gene ID</b>           | 6926                                               |
| <b>SwissProt ID</b>      | O15119.                                            |

## Application

|                         |                                            |
|-------------------------|--------------------------------------------|
| <b>Dilution Ratio</b>   | WB: 1:500-1:1000 IP: 1:20                  |
| <b>Molecular Weight</b> | Calculated MW: 79 kDa; Observed MW: 79 kDa |

## Background

---

**Product Name: T Box Protein 3 Rabbit Monoclonal Antibody**  
**Catalog #: AMRe02671**

---




Transcriptional repressor involved in developmental processes. Probably plays a role in limb pattern formation.

**Research Area**

Epigenetics and Nuclear Signaling

**Image Data**



Western blot analysis of Tbx3 in 3T3 lysates using T Box Protein 3 antibody.

**Note**

For research use only.