

---

**Product Name: TCTP Rabbit Monoclonal Antibody****Catalog #: AMRe02676**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,ICC/IF,IP
<b>Reactivity</b>	Human,Mouse,Rat
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	0.68mg/ml. The concentration of this product may be batch-dependent.
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	WB 1:500-1:1000,ICC/IF 1:50-1:200,IP 1:20-1:50
<b>Molecular Weight</b>	Calculated MW: 20 kDa; Observed MW: 20 kDa

**Antigen Information**

<b>Gene Name</b>	TPT1
<b>Alternative Names</b>	Translationally-controlled tumor protein; TCTP; Fortilin; Histamine-releasing factor; HRF; p23
<b>Gene ID</b>	7178
<b>SwissProt ID</b>	P13693
<b>Immunogen</b>	A synthetic peptide of human TCTP


**Background**

Involved in calcium binding and microtubule stabilization.


## Research Area

Immunology

## Image Data



Western blot analysis of TCTP in K562, rat Brain, C6, 3T3, HeLa lysates using TCTP antibody.



Immunocytochemistry analysis of TCTP(green) in HeLa using TCTP antibody, and DAPI(blue)