

**Product Name: Transmembrane Protein 161B Rabbit
Monoclonal antibody
Catalog #: AMRe02693**



Summary

Production Name	Transmembrane Protein 161B Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	TMEM161B
Alternative Names	Transmembrane protein 161B
Gene ID	153396
SwissProt ID	Q8NDZ6.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 55 kDa; Observed MW: 55 kDa

Background

**Product Name: Transmembrane Protein 161B Rabbit
Monoclonal antibody
Catalog #: AMRe02693**




Essential for maintaining normal cardiac rhythm in the developing heart and for neonatal survival (By similarity).

Research Area

Cell Biology

Image Data



Western blot analysis of TMEM161B in K562, HeLa lysates using Transmembrane Protein 161B antibody.

Note

For research use only.