

**Product Name: Thioredoxin Like Protein 1 Rabbit
Monoclonal Antibody
Catalog #: AMRe02731**

Summary

Production Name	Thioredoxin Like Protein 1 Rabbit Monoclonal Antibody
Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	TXNL1
Alternative Names	Txl; TXNL; TRP32; TXL-1; HEL-S-114
Gene ID	9352
SwissProt ID	O43396.

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
Molecular Weight	Calculated MW: 32 kDa; Observed MW: 32 kDa

**Product Name: Thioredoxin Like Protein 1 Rabbit
Monoclonal Antibody
Catalog #: AMRe02731**


Background

Active thioredoxin with a redox potential of about -250 mV.


Research Area

Signal Transduction


Image Data



Western blot analysis of TXNL1 in K562, rat Brain, C6, 3T3, HeLa lysates using Thioredoxin Like Protein 1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using TXNL1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



**Product Name: Thioredoxin Like Protein 1 Rabbit
Monoclonal Antibody
Catalog #: AMRe02731**



Immunocytochemistry analysis of TXNL1(green) in Hela using TXNL1 antibody,and DAPI(blue)

Note

For research use only.