

Product Name: UCHL5 Rabbit Monoclonal Antibody
Catalog #: AMRe02740



Summary

Production Name	UCHL5 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF, IP
Reactivity	Human, Mouse, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	UCHL5
Alternative Names	UCH37; CGI-70; INO80R; UCH-L5
Gene ID	51377
SwissProt ID	Q9Y5K5.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200 IP: 1:20
Molecular Weight	Calculated MW: 38 kDa; Observed MW: 38 kDa

Background

Product Name: UCHL5 Rabbit Monoclonal Antibody
Catalog #: AMRe02740




Protease that specifically cleaves 'Lys-48'-linked polyubiquitin chains.


Research Area

Cell Biology


Image Data



Western blot analysis of UCH37 in rat spleen lysates using UCHL5 antibody.



Western blot analysis of UCH37 in 3T3, Hela lysates using UCH37 antibody.



Immunocytochemistry analysis of UCH37(green) in Hela using UCH37 antibody, and DAPI(blue)

Product Name: UCHL5 Rabbit Monoclonal Antibody
Catalog #: AMRe02740



Note

For research use only.