

Product Name: HAUS7 Rabbit Monoclonal Antibody
Catalog #: AMRe02741



Summary

Production Name	HAUS7 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	HAUS7
Alternative Names	UIP1; UCHL5IP
Gene ID	55559
SwissProt ID	Q99871.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 40 kDa; Observed MW: 40 kDa

Background

Contributes to mitotic spindle assembly, maintenance of centrosome integrity and completion of cytokinesis as part of the

Product Name: HAUS7 Rabbit Monoclonal Antibody
Catalog #: AMRe02741




HAUS augmin-like complex.

Research Area

Cell Biology

Image Data



Western blot analysis of UCHL5IP in K562 lysates using HAUS7 antibody.

Note

For research use only.