

Product Name: USO1 Rabbit Monoclonal Antibody
Catalog #: AMRe02759



Summary

Production Name	USO1 Rabbit Monoclonal Antibody
Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	USO1
Alternative Names	TAP; VDP; P115
Gene ID	8615
SwissProt ID	O60763.

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:20-1:50
Molecular Weight	Calculated MW: 108 kDa; Observed MW: 108 kDa

Background

Product Name: USO1 Rabbit Monoclonal Antibody
Catalog #: AMRe02759




General vesicular transport factor required for intercisternal transport in the Golgi stack; it is required for transcytotic fusion and/or subsequent binding of the vesicles to the target membrane. May well act as a vesicular anchor by interacting with the target membrane and holding the vesicular and target membranes in proximity.

Research Area

Signal Transduction

Image Data



Western blot analysis of USO1 in 3T3 lysates using USO1 antibody.

Note

For research use only.