

Product Name: Villin Rabbit Monoclonal Antibody
Catalog #: AMRe02760



Summary

Production Name	Villin Rabbit Monoclonal Antibody
Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	VIL1
Alternative Names	VIL; VIL1; Villin1
Gene ID	7429
SwissProt ID	P09327.

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
Molecular Weight	Calculated MW: 93 kDa; Observed MW: 93 kDa

Background

Product Name: Villin Rabbit Monoclonal Antibody
Catalog #: AMRe02760




Ca(2+)-regulated actin-binding protein.

Research Area

Signal Transduction

Image Data



Western blot analysis of Villin in MCF-7 lysates using Villin antibody.

Note

For research use only.